

CRY ❄️ TEC 360°

MULTI-FUNCTIONAL CRYOLIPOLYSIS SLIMMING MACHINE




- ▶ **MORE COMFORTABLE**
- ▶ **MORE EFFECTIVE**
- ▶ **MORE TREATMENT AREAS**

FAT REDUCTION VS WEIGHT LOSS


It may be confused to hear that reducing fat is not the same as losing weight. It's both about what happens to fat cells. Losing weight makes fat cells smaller, and losing weight does not make fat cells go away. The number of fat cells in our bodies typically becomes a fixed number sometime in our teenage years. Losing or gaining weight generally does not increase or decrease the number of fat cells. Instead, it changes the size of fat cells.

FAT CELLS AFTER TYPICAL WEIGHT LOSS




After dieting and weight loss, fat cells get smaller but the number of fat cells stays the same. When weight is gained, fat cells get bigger.

FAT CELLS AFTER THE CRYOLIPOLYSIS PROCEDURE



The Cryolipolysis fat-freezing procedure actually reduces the number of fat cells in treated areas by about 20% to 25%. These fat cells are naturally processed and eliminated, so they won't migrate to another area of your body.

CRYOTEC360 VS OLD TECHNOLOGY



NEW GENERATION 360 DEGREE HANDPIECE

High-definition touch screen
Screen size 3.5inch

Medical soft silicone material



Replaceable contours



Change contours with ease

Contours are securely fastened to the handle quickly ,
very convenient for different areas treatment



10 different sizes handle cups, perfect for whole body
fat freezing treatment



Handle A 2 cups



Handle B 4 cups



235*65*47mm



255*65*55mm



180*65*50mm



190*65*55mm



210*65*70mm



220*65*80mm



Handle C 3 cups



125*50*45mm



145*50*70mm



135*50*50mm



Handle D

83*30*50mm

CRY ❄️ TEC 360°



Cavitation handle for arm,
leg and back fat break



RF handle for body shaping,skin
tightening after fat reduction



RF handle for double
chin treatment and face lifting





Easy to Operate

Intelligent design software, treatment parameters will be setted automatically according to the thickness of the fat.



WHY CRYOTEC 360 IS DIFFERENT FROM OTHER MACHINES



Shorter Treatment Time

Revolutionary cup design and four handpieces working simultaneous, save staff cost, save treatment time.



Target More Unwanted Fat Per Treatment

360 degree handpiece, 100% effective cooling area. Surrounding cooling action in entire applicator. Treatment result is much better.



Treatment More Stubborn Area


Interchangeable applicators with 10 different size contours allows us to treat small and large areas of fat on the neck, arm, body and thighs.



100% Safety Treatment

Professional temperature sensors are installed inside the handle to monitor the cooling temperature in real time to ensure the safety of treatment.

UNIQUE TEC COOLING SYSTEM



The instrument is internally loaded with a TEC refrigeration system and an intelligent temperature monitoring sensor to ensure that the operating temperature of the machine is constant at 20 degrees Celsius. Make sure that the water temperature remains the same for 24 hours of continuous operation.

SUMMARY OF YOUR TREATMENT



Procedure Time

40-60 minutes per application estimated



Full recovery

Immediate



Number of Treatments

1-3



Back To Work

Immediate estimated



Sensitivity Period

Up to 7 days estimated



Risks & Complications

Slight numbness for up to 7 days



Anaesthetic

None



Results

Visible after 2-3 weeks



Duration Of Results

Long Term Results

BEFORE & AFTER



Before



After 6 weeks



Before



After 8 weeks



Before



After 4 weeks



Before



After 8 weeks



Before



After 9 weeks



Before



After 6 weeks